

## ► FMT Disposable Pulseoxymeter (SpO2) Sensors Specifications

FMT disposable pulse oximetry sensors are ideal for extended monitoring or minimizing the risk of infection where the risk of cross-contamination is high.

FMT disposable sensors are made of flexible shielded cable, high accuracy electronic components and 3M Microfoam™ or Medaplast™ fabric adhesive materials. These adhesive materials with medium adhesion maximize patient comfort and allow the sensor to be easily moved and quickly reapplied when repositioning.

FMT also offers velcro non-adhesive single patient use sensors. Due to the Velcro design, the repositioning of the sensor is possible for many times without destroying the sensor body. This sensor is very helpful to users especially during long term monitorization.

FMT disposable sensors have special design that ensures easy and accurate sensor placement and functionality. All FMT disposable sensors are Latex-free and meet industry standards for biocompatibility.

FMT offers four sizes of comfortable and accurate disposable sensors designed to fit adult, pediatric, infant and neonatal patients.


FMT disposable sensors are compatible with major leading manufacturers of patient monitoring and pulse oximetry systems.

FMT disposable pulse oximetry sensors have high signal accuracy. Accuracy is  $\pm 1\%$  between 90%–100% and  $\pm 2\%$  between 70%–89%. Added Faraday shielding feature reduces the effect of magnetic interference from the environment and increases the accuracy of the measurement.

FMT disposable pulse oximetry sensors are supplied sterile or non-sterile and sold in 24 pieces per pack.

FMT disposable pulse oximetry sensors are compatible with PWTT technology.

## ► Disposable SpO2 Sensors with Microfoam type material








# FMT

## DISPOSABLE PULSE OXIMETRY (SPO2) SENSORS

### ► Disposable SpO2 Sensors with Fabric type material



Sensor Type	Patient Type	Patient Weight	Cable Length	Preferred Application Area
FMT-DAF/XXX-Y	Adult	>30 kg (>66 lbs)	45 cm	Index finger
FMT-DPF/XXX-Y	Pediatric	10-50 kg (22-110 lbs)	45 cm	Index finger
FMT-DIF /XXX-Y	Infant	1-20 kg (2.2-44 lbs)	90 cm	Big Toe
FMT-DNF/XXX-Y	Neonatal/Adult	<3kg or >30 kg (<6.6 lbs or > 66lbs)	90 cm	Lateral aspect of foot or Index finger (for adult)

Compatible Brand	Patient Type			
	Adult	Pediatric	Infant	Neonatal / Adult
BCI*	FMT-DAF / BCI - Y	FMT-DPF / BCI - Y	FMT-DIF / BCI - Y	FMT-DNF / BCI - Y
CSI*	FMT-DAF / CSI - Y	FMT-DPF / CSI - Y	FMT-DIF / CSI - Y	FMT-DNF / CSI - Y
DATEX*	FMT-DAF / DTX - Y	FMT-DPF / DTX - Y	FMT-DIF / DTX - Y	FMT-DNF / DTX - Y
GE OXY TIP (UN Connector)*	FMT-DAF / OXY - Y	FMT-DPF / OXY - Y	FMT-DIF / OXY - Y	FMT-DNF / OXY - Y
MEK*	FMT-DAF / MEK - Y	FMT-DPF / MEK - Y	FMT-DIF / MEK - Y	FMT-DNF / MEK - Y
NELLCOR	FMT-DAF / NLC - Y	FMT-DPF / NLC - Y	FMT-DIF / NLC - Y	FMT-DNF / NLC - Y
NELLCOR OXIMAX	FMT-DAF / NLO - Y	FMT-DPF / NLO - Y	FMT-DIF / NLO - Y	FMT-DNF / NLO - Y
NIHON KOHDEN	FMT-DAF / NKN - Y	FMT-DPF / NKN - Y	FMT-DIF / NKN - Y	FMT-DNF / NKN - Y
NONIN*	FMT-DAF / NON - Y	FMT-DPF / NON - Y	FMT-DIF / NON - Y	FMT-DNF / NON - Y
NOVAMETRIX*	FMT-DAF / NVM - Y	FMT-DPF / NVM - Y	FMT-DIF / NVM - Y	FMT-DNF / NVM - Y

\* Manufactured upon order

Y: "Y" is for wrap material. Please add "M" for Microfoam, "F" for fabric material.

Example: FMT-DAF/NLC-F means fabric material, adult size, disposable SpO2 sensor for Nellcor.





# FMT

## DISPOSABLE PULSE OXIMETRY (SPO2) SENSORS

### ► FMT-DVF Series single patient use SpO2 sensors

FMT single patient use sensor	Patient Type	Patient Weight	Cable Length	Preferred Application Area
FMT-DVF/XXX	Neonatal/Adult	<3kg or >30 kg (<6.6 lbs veyá >66lbs)	90 cm	Lateral aspect of foot or Index finger (for adult)



FMT-DVF/NLC



FMT-DVF/OXY

Compatible Brand	Patient Type
	Neonatal / Adult
BCI*	FMT-DVF / BCI
CSI*	FMT-DVF / CSI
DATEX*	FMT-DVF / DTX
GE OXY TIP (UN Connector)*	FMT-DVF / OXY
MEK*	FMT-DVF / MEK
NELLCOR	FMT-DVF / NLC
NELLCOR OXIMAX	FMT-DVF / NLO
NIHON KOHDEN	FMT-DVF / NKN
NONIN*	FMT-DVF / NON
NOVAMETRIX*	FMT-DVF / NVM

\* Manufactured upon order



EN ISO 9001:2008  
EN ISO 13485:2003

FMT® is registered trademark of Metko Ltd. Other products, logos and company names mentioned herein maybe trademarks of their respective owners and are not to be taken as an endorsement or affiliation with FMT®. The information in this instruction insert has been carefully checked and it is believed to be accurate. In the interest of continued product development, METKO reserves the right to make changes and improvements to this insert and the product it described any time, without notice or obligation.



#### HEADQUARTER & FACTORY

- İvedik O.S.B. Ağaç İşleri Sanayi Sit.  
1354. Cad. 1358.Sok. No:9  
06378 Yenimahalle

ANKARA / TÜRKİYE

Tel (pbx) : +90 312 387 12 46  
Fax : +90 312 387 12 51

#### BRANCH

- Atatürk Mah. Sedef Cad. 38 Ada  
Ata Plaza 3-5 Blokları Kat:2 No:57  
34758 Ataşehir

İSTANBUL / TÜRKİYE

Tel (pbx) : +90 216 455 91 21  
Fax : +90 216 455 48 24

#### INTERNET

- <http://www.metkoltd.com>  
[metko@metkoltd.com](mailto:metko@metkoltd.com)