



LISCLARE™
Freedom to Care

Contact

For more information and to request a free sample please contact:

📞 GB: 0151 495 5750 NI: 028 9079 4000 ROI: 01 633 2020

✉️ customerservice@lisclare.com

PRACTICAL benefits

A high safe working load of up to 500kg caters for most patients



The wide range of styles covers most requirements including high-back, amputee, hygiene and sit-to-stand.



The versatility of the MySling range means that there is no need to look elsewhere for costly, specialist slings.

Can be used in hygiene environment, hydrotherapy areas, showers, baths and pools

Available with clip or loop attachments



SAFETY

COMFORT

VALUE

MySling®
Lasting Disposable Slings

1 SAFETY in the real world


Laboratory testing does not always guarantee safety in the real world, which is why our disposable slings feature a threaded **RipStop** mesh to prevent tearing during use.

Everyday events during sling usage can create a weak point in the fabric. These events include accidental snagging on furniture or simply applying excess force to the fabric by pulling or gripping the sling when manouvering the client. Moisture caused by incontinence or hygiene activities can also degrade most non-woven fabrics.



Our **RipStop** fabric prevents tearing or splitting caused by:

- Snagging
- Wrenching
- Moisture



2 No compromise on COMFORT



Our multi-layer, stitch-bonded fabric is soft, breathable and comfortable. MySling fabric is designed to distribute the patient's weight evenly and eliminate pressure points.

Moving to MySling allows you to improve infection control without reducing patient comfort.

- ✓ Soft
- ✓ Breathable
- ✓ Even weight distribution

3 LOWER Monthly Costs

Slings which are safe to reuse lead to substantial monthly savings:

"After switching from conventional disposable clip slings to MySling, our annual spend reduced from £90k to £52k per year, an annual saving of 42%."

Procurement Manager,
NHS Trust Hospital, North Midlands

Reuse with
confidence

The unique safety features of the MySling range mean that the same sling can be reused for up to 6 months. The danger of reusing conventional disposables is that, if they become moist, or have been wrenched or snagged, they may not be safe for the patient to use a second time.

RipStop technology ensures that a MySling product is safe to use over and over again, resulting in reduced monthly usage and therefore real monthly cost savings.



Reduce waste