

Ultrasound Anytime,  
**Anywhere.**  
**ACUSON P500**





Enabling  
ultrasound  
imaging anytime,  
anywhere.

# Ultrasound Anytime, Anywhere

Siemens Healthineers engineered the compact and powerful ACUSON P500™, a portable ultrasound system designed to equip clinicians with a powerful lineup of clinical applications right at the bedside.

Focused entirely on the user, the ACUSON P500 is your go-to platform to evaluate clinical manifestations, facilitate accurate diagnoses, and to assist procedures in emergency and critical care.

The ACUSON P500, with a smaller footprint in a crowded technology environment, also offers an IntraCardiac Echocardiography (ICE) Edition

that integrates the imaging capabilities of the ACUSON AcuNav™ family of catheters to the ACUSON P500 system. The utilization of the 10F & 8F ACUSON AcuNav ultrasound catheters provide real-time visualization of cardiac anatomy, catheters and devices within the heart and creates 3D contours of the heart chambers with the SoundStar® catheter. In addition, the bedside remote control capability allows the clinician quick access to image optimization.

With the ACUSON P500, Siemens Healthineers is enabling quality care with ultrasound imaging anytime, anywhere.



## See better

Overcome the challenges of motion artifacts and achieve high quality images.



## Scan faster

Ability to greatly improve workflow powered by an innovative and user friendly interface with integrated touch technology and quick boot-up times.




## Go further

Offer exceptional detail and quality of information with over 13 innovative transducers, advanced intra-cardiac echo imaging capabilities, needle guidance solutions, and much more.

# Advanced Technology to See Better, Scan Faster & Go Further

The ACUSON P500 enables reliable and fast imaging by offering new and pioneering high-end imaging technology. It provides standardized, high-quality images with more depth, clearer border definition and more tissue variation.




Scan rapidly through lower extremity veins with excellent image clarity and reduced image noise with Dynamic Persistence.



Dynamic Persistence and Auto Flash Artifact Suppression reduce tissue motion and color Doppler artifact, such as the lower extremity vessels.



Greatly reduce flash artifact in highly pulsatile vessels, such as the hepatic veins with Auto Flash Artifact Suppression.




Touch interface allows users to more easily interact with the imaging parameters on the screen for increased workflow efficiency.



Auto Flash Artifact Suppression reduces color Doppler flash while simultaneously enhancing color sensitivity in the carotid artery.




Motion correction technology improves diagnostic confidence while performing dynamic scans of Parasternal Long-Axis (PLAX) view.




Improve tissue differentiation and border detection, such as in this breast cyst, with Dynamic TCE™ tissue contrast enhancement technology.




Excellent contrast resolution and visualization of subtle tissue differences, such as in the kidney, with Advanced SieClear™ spatial compounding.




Anticipate and execute your exam based on customizable programs using eSieScan™ workflow protocols such as in this annotation and measurement of Apical 3 Chamber with color.



Auto Flash Artifact Suppression allows you to scan rapidly without flash artifact, such as in this color Doppler of a thyroid nodule.



Motion correction technology improves diagnostic confidence while performing dynamic scans of joints in motion, such as this shoulder.




Configure variety of stress test protocols including this Apical 2 chambers stress test using the ACUSON P500.

## See Better

The ACUSON P500 improves image quality with exclusive motion correction and a diverse range of applications and features designed specifically for your area of expertise. With Siemens' innovative technology, image quality is paramount.

Therapeutic exams  
less complicated  
and more  
standardized



Reduce motion artifacts and noise while simultaneously enhancing color sensitivity with Dynamic Persistence and Auto Flash Artifact Suppression technology.



Visualize subtle tissue variations with Dynamic TCE™ tissue contrast enhancement technology, enhance anatomical detail with Advanced SieClear™ spatial compounding and optimize images with one button, employing TEQ™ ultrasound technology.



Leverage Clarify™ vascular enhancement (VE) technology that uniquely utilizes power Doppler flow information to enhance B-mode imaging.






## Scan Faster

The ACUSON P500 offers a wide range of innovative applications that make diagnostic and therapeutic exams less complicated and more standardized. Exclusive automated solutions available with the ACUSON P500 improve scan time, expedite workflow and increase operational efficiency.

Perform a wide range of exams with the Stress Echo application, which creates user-defined protocols and offers flexible review formats and visual analysis.

Improve accuracy and reproducibility with the eSie Measure™ workflow acceleration package, the first application that provides automated measurements for routine echo exams.

Reduce training time, variability in skill and operation errors with eSieScan™ Protocol, a unique semi-automated guided workflow designed to automate and simplify routine tasks.





Work more efficiently,  
consistently and  
clearly with powerful  
applications.



## Go Further

The ACUSON P500, with its unmatched diagnostic and workflow solutions, enables you to work more efficiently, consistently and clearly with powerful applications that enable you to increase patient throughput.

Improve visibility with Needle Visualization Enhancement that employs unique multi-beam spatial compounding and other optimization techniques while still maintaining excellent image quality, especially at steeper angles.

Utilize a wide range of anatomy specific transducers that cover multiple applications.

Reduce noise for image clarity and minimize artifact caused by ablation technique using the newly engineered SwiftLink™ catheter connector for ACUSON P500.




## Go Further

Your reputation as a healthcare provider and your financial success rely on best-in-class staff expertise, exceptional system performance, and seamless clinical workflows. Designed to adapt to highest standards of care – from technical compliance standards to system usability – with the ACUSON P500 you can rest assured that the equipment you invest in today will be up-to-date and ready to integrate the most advanced technologies for years to come. This allows you to manage costs without compromising the quality of care.

### eSieLink Remote Assistance Technology

A new service tool for your ultrasound systems, eSieLink supports your clinical staff with fast answers to technical and clinical questions. Secure remote desktop sharing allows them to communicate in real-time with technical experts from Siemens Healthineers to resolve issues right away and to receive additional training support. eSieLink remote assistance technology helps to improve system uptime and promotes seamless clinical workflows.



# Feature Overview

## ACUSON P500 Ultrasound System

The portable ultrasound system is designed to provide clinicians with an innovative and diverse range of applications and features at the point-of-care so they are able to see better, scan faster and go further.

## ACUSON P500 Ultrasound System, ICE Edition

The Intracardiac Echocardiography Edition integrates the imaging capabilities of Siemens AcuNav™ ultrasound catheter technologies with the ACUSON P500 ultrasound system to provide real-time visualization of cardiac anatomy.

Options	ACUSON P500™ ultrasound system	ACUSON P500™ ultrasound system, ICE Edition
Motion Correction Technology	Standard	Standard
Tissue Harmonic Imaging (THI)	Standard	Standard
DTI™ Doppler tissue imaging capability	Standard	Standard
Advanced SieClear™ Spatial Compounding	Standard	Standard
Wireless Data Transfer	Standard	Standard
Dynamic TCE™ tissue contrast enhancement technology	Standard	Standard
Clarify™ vascular enhancement (VE) technology	Standard	Standard
Cardiac M&R	Standard	Standard
Stress Echo Imaging	Optional	Optional
SieScape™ panoramic imaging	Optional	Optional
syngo® Arterial Health Package (AHP)	Optional	Optional
Pencil Transducers	Optional	Optional
eSie Measure™ workflow acceleration package	Optional	Optional
Enhanced Needle Visualization	Optional	Optional
eSieScan™ workflow protocols	Optional	Optional
Barcode Reader	Optional	Optional
DICOM 3.0 Connectivity (Option: DICOM Structured Reporting)	Optional	Optional
Intracardiac Echocardiography	N/A	Standard
SwiftLink™ catheter connector	N/A	Optional
JoyStick	N/A	Optional
CartoSound® Communication	N/A	Optional

The products/features mentioned in this document may not be commercially available in all countries. Due to regulatory reasons their future availability cannot be guaranteed. Please contact your local Siemens organization for further details. Standalone clinical images may have been cropped to better visualize pathology.

ACUSON P500, AcuNav, Clarify, DTI, eSie Measure, eSie Scan, SieClear, SieScape, SwiftLink, TCE and TEQ are trademarks of Siemens Medical Solutions USA, Inc.

CartoSound and SoundStar are registered trademarks of Biosense Webster, part of the Johnson & Johnson Family of Companies.

*syngo*® is a registered trademark owned by Siemens Healthcare GmbH.

*To learn more about ACUSON P500 Ultrasound Systems, please visit us at [siemens.com/ultrasound](https://www.siemens.com/ultrasound).*

#### **Siemens Healthineers Headquarters**

Siemens Healthcare GmbH  
Henkestr. 127  
91052 Erlangen  
Germany  
Phone: +49-9131-84-0  
[siemens.com/healthineers](https://www.siemens.com/healthineers)

#### **Legal Manufacturer**

Siemens Medical Solutions USA, Inc. Ultrasound  
685 E. Middlefield Road  
Mountain View, CA 94043 USA  
Phone: 1-888-826-9702  
[siemens.com/healthcare](https://www.siemens.com/healthcare)